

Product datasheet for **MR216938A1V**

Mouse Cxcl3 (NM_203320) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: Cxcl3

Synonyms: Dci; Dcpl; Gm1960

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR216938 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCCTCCACCTGCCGGCTCCTCAGTGCTGCACTGGTCCTGCTGCTGCTGGCCACCAACCACC
AGGTACAGGGGCTGTTGTGGCCAGTGAGCTGCGCTGTCAGTGCCTGAACACCTACCAAGGTTGATT
TGAGACCATCCAGAGCTTGACGGTGACGCCCCAGGACCCCACTGCACCCAGACAGAAGTCATAGCCACT
CTCAAGGATGGTCAAGAAGTTGCCTCAACCCCAAGGCCCAAGGCTTCAGATAATCATCAAGAAGATAC
TGAAGAGCGGCAAGTCCAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216938 protein sequence
Red=Cloning site Green=Tags(s)

MAPPTCRLLSAALVLLLLLATNHQATGAVVASELRQCCLNTLPRVDFETIQSLTVTPPGPHCTQTEVIAT
LKDQGEVCLNPQGPRLQIIKKILKSGKSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse

Serotype: AAV-2

ACCN: NM_203320

ORF Size: 300 bp



Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_203320.2</u>
RefSeq Size:	1018 bp
RefSeq ORF:	303 bp
Locus ID:	330122
UniProt ID:	<u>Q6W5C0</u>
Cytogenetics:	5 E1
MW:	10.7 kDa