

Product datasheet for **MR216958L3V**

Cd302 (NM_025422) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd302 (NM_025422) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd302
Synonyms:	1110055L24Rik; A1159627
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025422
ORF Size:	657 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216958).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025422.3 , NP_079698.2
RefSeq Size:	1257 bp
RefSeq ORF:	660 bp
Locus ID:	66205
UniProt ID:	Q9DCG2
Cytogenetics:	2 C1.1
Gene Summary:	Potential multifunctional C-type lectin receptor that may play roles in endocytosis and phagocytosis as well as in cell adhesion and migration.[UniProtKB/Swiss-Prot Function]



[View online »](#)