

Product datasheet for **MR216965L4V**

Rhd (NM_011270) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rhd (NM_011270) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rhd
Synonyms:	Rh; Rhced; Rh11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011270
ORF Size:	1254 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216965).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011270.3 , NP_035400.3
RefSeq Size:	1530 bp
RefSeq ORF:	1257 bp
Locus ID:	19746
UniProt ID:	Q8CF94
Cytogenetics:	4 67.13 cM
Gene Summary:	May be part of an oligomeric complex which is likely to have a transport or channel function in the erythrocyte membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)