

Product datasheet for **MR216986L3V**

Tmx2 (NM_025868) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmx2 (NM_025868) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmx2
Synonyms:	2310042M24Rik; AA589631; Txndc14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025868
ORF Size:	885 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216986).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025868.4 , NP_080144.1
RefSeq Size:	2287 bp
RefSeq ORF:	888 bp
Locus ID:	66958
UniProt ID:	Q9D710
Cytogenetics:	2 D



[View online »](#)

Gene Summary:

Endoplasmic reticulum and mitochondria-associated protein that probably functions as a regulator of cellular redox state and thereby regulates protein post-translational modification, protein folding and mitochondrial activity. Indirectly regulates neuronal proliferation, migration, and organization in the developing brain.[UniProtKB/Swiss-Prot Function]