

Product datasheet for **MR216999L3V**

Arhgap28 (NM_172964) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgap28 (NM_172964) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgap28
Synonyms:	AU044757; AW550892; E130310N06
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172964
ORF Size:	2187 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216999).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172964.4 , NP_766552.3
RefSeq Size:	5317 bp
RefSeq ORF:	2190 bp
Locus ID:	268970
UniProt ID:	Q8BN58
Cytogenetics:	17 E1.1-E1.2
Gene Summary:	GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.[UniProtKB/Swiss-Prot Function]



[View online »](#)