

Product datasheet for **MR217018L3V**

Aen (NM_001162939) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aen (NM_001162939) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aen
Synonyms:	2700083B06Rik; Isg2011
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001162939
ORF Size:	1008 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217018).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001162939.1
RefSeq Size:	2466 bp
RefSeq ORF:	897 bp
Locus ID:	68048
UniProt ID:	Q9CZ19
Cytogenetics:	7 D2



[View online »](#)

Gene Summary:

Exonuclease with activity against single- and double-stranded DNA and RNA. Mediates p53-induced apoptosis. When induced by p53 following DNA damage, digests double-stranded DNA to form single-stranded DNA and amplifies DNA damage signals, leading to enhancement of apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]