

Product datasheet for **MR217040L4V**

Prmt3 (NM_133740) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prmt3 (NM_133740) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Prmt3
Synonyms:	2010005E20Rik; 2410018A17Rik; AL033309; Hrmt1I3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133740
ORF Size:	1584 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217040).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133740.2 , NP_598501.1
RefSeq Size:	2495 bp
RefSeq ORF:	1587 bp
Locus ID:	71974
UniProt ID:	Q922H1
Cytogenetics:	7 B4
Gene Summary:	Methylates (mono and asymmetric dimethylation) the guanidino nitrogens of arginyl residues in some proteins.[UniProtKB/Swiss-Prot Function]



[View online »](#)