

Product datasheet for **MR217042L4V**

Ccdc39 (NM_026222) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccdc39 (NM_026222) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ccdc39
Synonyms:	4921507O14Rik; AI447426; AI844750; b2b1304Clo; b2b1735Clo; b2b2025.1Clo; D3Ertd789e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026222
ORF Size:	2814 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217042).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026222.2 , NP_080498.1
RefSeq Size:	3616 bp
RefSeq ORF:	2814 bp
Locus ID:	51938
UniProt ID:	Q9D5Y1
Cytogenetics:	3 A3



[View online »](#)

Gene Summary:

Required for assembly of dynein regulatory complex (DRC) and inner dynein arm (IDA) complexes, which are responsible for ciliary beat regulation, thereby playing a central role in motility in cilia and flagella. Probably acts together with CCDC40 to form a molecular ruler that determines the 96 nanometer (nm) repeat length and arrangements of components in cilia and flagella. Not required for outer dynein arm complexes assembly.[UniProtKB/Swiss-Prot Function]