

Product datasheet for **MR217046L3V**

Epb4111 (NM_013510) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Epb4111 (NM_013510) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Epb4111
Synonyms:	4.1N; Epb4.111; mKIAA0338; NBL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013510
ORF Size:	2637 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217046).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013510.3
RefSeq Size:	6298 bp
RefSeq ORF:	2640 bp
Locus ID:	13821
UniProt ID:	Q9Z2H5
Cytogenetics:	2 77.39 cM
Gene Summary:	May function to confer stability and plasticity to neuronal membrane via multiple interactions, including the spectrin-actin-based cytoskeleton, integral membrane channels and membrane-associated guanylate kinases.[UniProtKB/Swiss-Prot Function]



[View online »](#)