

Product datasheet for **MR217047L4V**

Pdpr (NM_198308) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Pdpr
Synonyms:	4930402E16Rik; AU018056; mKIAA1990
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198308
ORF Size:	2634 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR217047).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_198308.1](#), [NP_938050.1](#)

RefSeq Size: 2904 bp

RefSeq ORF: 2637 bp

Locus ID: 319518

UniProt ID: [Q7TSQ8](#)

Cytogenetics: 8 E1



Gene Summary:

Decreases the sensitivity of PDP1 to magnesium ions, and this inhibition is reversed by the polyamine spermine.[UniProtKB/Swiss-Prot Function]