

Product datasheet for **MR217054L3V**

Tcn2 (NM_001130458) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tcn2 (NM_001130458) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tcn2
Synonyms:	AW208754; Tcn-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130458
ORF Size:	1290 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217054).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130458.1 , NP_001123930.1
RefSeq Size:	2058 bp
RefSeq ORF:	1293 bp
Locus ID:	21452
UniProt ID:	O88968
Cytogenetics:	11 2.76 cM
Gene Summary:	Primary vitamin B12-binding and transport protein. Delivers cobalamin to cells. [UniProtKB/Swiss-Prot Function]



[View online »](#)