

## Product datasheet for MR217065L3V

### Sec61g (NM\_001109971) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sec61g (NM_001109971) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sec61g
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001109971
ORF Size:	204 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217065).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001109971.1</a> , <a href="#">NP_001103441.1</a>
RefSeq Size:	779 bp
RefSeq ORF:	207 bp
Locus ID:	20335
UniProt ID:	<a href="#">P60060</a>
Cytogenetics:	11 A2
Gene Summary:	Component of SEC61 channel-forming translocon complex that mediates transport of signal peptide-containing precursor polypeptides across endoplasmic reticulum (ER). [UniProtKB/Swiss-Prot Function]



[View online »](#)