

## Product datasheet for **MR217073L3V**

### H13 (NM\_010376) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | H13 (NM_010376) Mouse Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | H13  |
| Synonyms:                 | 1200006O09Rik; 4930443L17Rik; 5031424B04Rik; AV020344; H-13; Hm13; PSL3; Spp   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_010376  |
| ORF Size:                 | 1134 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR217073).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_010376.4</a> , <a href="#">NP_034506.1</a>  |
| RefSeq Size:              | 5202 bp  |
| RefSeq ORF:               | 1137 bp  |
| Locus ID:                 | 14950  |
| UniProt ID:               | <a href="#">Q9D8V0</a>   |
| Cytogenetics:             | 2 75.41 cM   |



[View online »](#)

**Gene Summary:**

Catalyzes intramembrane proteolysis of some signal peptides after they have been cleaved from a preprotein, resulting in the release of the fragment from the ER membrane into the cytoplasm. Required to generate lymphocyte cell surface (HLA-E) epitopes derived from MHC class I signal peptides. Involved in the intramembrane cleavage of the integral membrane protein PSEN1. Cleaves the integral membrane protein XBP1 isoform 1 in a DERL1/RNF139-dependent manner (By similarity). May play a role in graft rejection (PubMed:9354467). [UniProtKB/Swiss-Prot Function]