

Product datasheet for **MR217079L4V**

Cebpe (NM_207131) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cebpe (NM_207131) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cebpe
Synonyms:	C/EBPe; CRP1; Gm294
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207131
ORF Size:	843 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217079).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207131.1 , NP_997014.1
RefSeq Size:	1238 bp
RefSeq ORF:	846 bp
Locus ID:	110794
UniProt ID:	Q6PZD9
Cytogenetics:	14 27.9 cM


[View online »](#)

Gene Summary:

Transcriptional activator. C/EBP are DNA-binding proteins that recognize two different motifs: the CCAAT homology common to many promoters and the enhanced core homology common to many enhancers. Required for the promyelocyte-myelocyte transition in myeloid differentiation.[UniProtKB/Swiss-Prot Function]