

Product datasheet for **MR217106L3V**

Hnrnpf (NM_001166429) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrnpf (NM_001166429) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrnpf
Synonyms:	4833420I20Rik; AA407306; Hnrpf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001166429
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217106).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001166429.1 , NP_001159901.1
RefSeq Size:	2197 bp
RefSeq ORF:	1248 bp
Locus ID:	98758
UniProt ID:	Q9Z2X1
Cytogenetics:	6 F1



[View online »](#)

Gene Summary:

Component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complexes which provide the substrate for the processing events that pre-mRNAs undergo before becoming functional, translatable mRNAs in the cytoplasm. Plays a role in the regulation of alternative splicing events. Binds G-rich sequences in pre-mRNAs and keeps target RNA in an unfolded state (By similarity).[UniProtKB/Swiss-Prot Function]