

## Product datasheet for **MR217108L3V**

### Strn4 (NM\_001039878) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Strn4 (NM_001039878) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Strn4
Synonyms:	C80611; ZIN; zinedin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001039878
ORF Size:	2259 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217108).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001039878.2</a> , <a href="#">NP_001034967.1</a>
RefSeq Size:	3541 bp
RefSeq ORF:	2262 bp
Locus ID:	97387
UniProt ID:	<a href="#">P58404</a>
Cytogenetics:	7 A2
Gene Summary:	Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.[UniProtKB/Swiss-Prot Function]



[View online »](#)