

## Product datasheet for **MR217119L3V**

### **Srgap3 (NM\_080448) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Srgap3 (NM_080448) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Srgap3
Synonyms:	AI452337; Arhgap14; D130026O08Rik; Gbi; MEGAP; mKIAA0411; WRP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080448
ORF Size:	3297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217119).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_080448.4</a> , <a href="#">NP_536696.4</a>
RefSeq Size:	8887 bp
RefSeq ORF:	3300 bp
Locus ID:	259302
UniProt ID:	<a href="#">Q812A2</a>
Cytogenetics:	6 E3
Gene Summary:	GTPase-activating protein for RAC1 and perhaps CDC42, but not for RhoA small GTPase. May attenuate RAC1 signaling in neurons (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)