

## Product datasheet for MR217139L3V

### Tmc3 (NM\_177695) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmc3 (NM_177695) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmc3
Synonyms:	9030203A06
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177695
ORF Size:	3390 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217139).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_177695.4</a> , <a href="#">NP_808363.3</a>
RefSeq Size:	4316 bp
RefSeq ORF:	3393 bp
Locus ID:	233424
UniProt ID:	<a href="#">Q7TQ69</a>
Cytogenetics:	7 D3
Gene Summary:	Probable ion channel.[UniProtKB/Swiss-Prot Function]



[View online »](#)