

## Product datasheet for MR217190

### Prkacb (NM\_001164200) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Prkacb (NM\_001164200) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Prkacb  
**Synonyms:** CbPKA; Pkacb  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR217190 representing NM\_001164200  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGGCAGCTCACAAAGAACTGTCCAGTGGCCAGCACAGTGGCAGCCCACCGCCCTCCAGAAGCTGGAGG  
 GGTTTCGCTAGCCGCCTGTTTCACAGGCACTCCAGAGCCACCGCACAAAGAGCACAGGGCAGCTCTGGAAGA  
 TGACGGCCTTCGGGCCTCTGAGCACACTGCCTCATGGGACAAATCAATGAAAGAGTTTCTAGCCAAAGCC  
 AAAGAAGACTTTCTGAGGAAATGGGAGAACCCTCCCCGAGTAATGCTGGCTTGAGGATTTTGGAGGA  
 AGAAAACCCTCGGGACGGGTTCCCTTGGAAAGATCATGTTGGTGAAGCATAAAGCCACTGAGCAGTACTA  
 CGCATGAAGATCTTAGACAAGCAGAAGTTGTTAAGCTGAAGCAAATAGAGCACACTCTGAATGAGAAG  
 AGAATCCTGCAGGCCGTGGAGTTCCCGTTCCTTGTGCGGCTGGAGTACTCTTTAAGGATAATTCTAATT  
 TATACATGGTTATGGAATACGTCCTCGGGGAGAGATGTTCTCACATCTGAGAAGAATTGGAAGGTTTCAG  
 TGAGCCCCACGCCGTTTCTATGCAGCCAGATTGTGCTAACATTTGAGTACCTTCATCCCTCGACCTC  
 ATCTACAGAGATCTCAAGCCGAAAACCTCTTAATTGACCACCAGGGTTACATCCAGGTCACAGATTTTCG  
 GGTTCCGCAAAAGAGTCAAGGGCAGGACATGGACATTGTGTGGCACCCAGAGTACCTGGCCCCGGAGAT  
 CATCCTCAGCAAGGTTACAATAAGCGGTGGACTGGTGGGCACTGGGCGTGTGATCTATGAGATGGCT  
 GCTGGTACCCCTCATTCTTTGCTGACCAGCCAATTCAGATCTATGAGAAGATTGCTCTGGAAAGGTTCC  
 GGTTCCCATCACACTTCAGCTCCGATCTCAAGGACCTTCTGCGGAACCTGCTGCAGGTGGATCTGACAAA  
 GCGATTCCGGAACCTGAAGAACGGCGTGTGATGACATAAAGACCCACAAGTGGTTTCCACAACTGACTGG  
 ATTGCTATTTATCAGAGAAAGGTTGAGGCTCCATTATACCAAAGTTTCAGAGGCTCTGGCGATACCAGCA  
 ACTTCGATGACTATGAAGAAGAAGAAATCCGTGTGTCTATAACAGAAAAATGTGGAAGGAATTTGTGA  
 ATTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR217190 representing NM\_001164200  
Red=Cloning site Green=Tags(s)

MAAHKELSSGQHSPTALQKLEGFASRLFHRHSRGTAEHRAALEDGLRASEHTASWDKSMKEFLAKA  
 KEDFLRKWENPPPSNAGLEDFERKKTGTGSFGRVMLVKHKATEQYYAMKILDKQKVVKLKQIEHTLNEK  
 RILQAVEFPFLVRLEYSFKDNSNLYMVMYVPGGEMFSLHRRIGRFSEPHARFYAAQIVLTFEYLHSLDL  
 IYRDLKPENLLIDHQGYIQVDFGFAKRVKGRWTLCGTPEYLAPEIILSKGYNKAVDWWALGVL IYEMA  
 AGYPPFFADQPIQIYEKIVSGKVRFP SHFSSDLKDLLRNLLQVDLTKRFGNLKNGVSDIKTHKWFATTDW  
 IAIYQRKVEAPFIPKFRGSGDTSNFDDYEEEEIRVSITEKCGKEFCF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9010\\_d03.zip](https://cdn.origene.com/chromatograms/mm9010_d03.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001164200


**ORF Size:** 1194 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001164200.1</a></u> , <u><a href="#">NP_001157672.1</a></u>
<b>RefSeq Size:</b>	4302 bp
<b>RefSeq ORF:</b>	1197 bp
<b>Locus ID:</b>	18749
<b>UniProt ID:</b>	<u><a href="#">P68181</a></u>
<b>Cytogenetics:</b>	3 H2
<b>MW:</b>	46.5 kDa
<b>Gene Summary:</b>	Mediates cAMP-dependent signaling triggered by receptor binding to GPCRs. PKA activation regulates diverse cellular processes such as cell proliferation, the cell cycle, differentiation and regulation of microtubule dynamics, chromatin condensation and decondensation, nuclear envelope disassembly and reassembly, as well as regulation of intracellular transport mechanisms and ion flux (PubMed:9368018). Regulates the abundance of compartmentalized pools of its regulatory subunits through phosphorylation of PJA2 which binds and ubiquitinates these subunits, leading to their subsequent proteolysis. Phosphorylates GPKOW which regulates its ability to bind RNA (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR217190