

## Product datasheet for **MR217206**

### **King1 (NM\_001102411) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	King1 (NM_001102411) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	King1
Synonyms:	King
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>MR217206 representing NM\_001102411  
Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGAAGCTCATTACTACTGCTCTCTGCTCCGACTCTGCTGACTTTAACACAGGGAGAAGAAGCGC  
AGGAAATTGACTGCAATGATGAGGCTGTATTTACGGCTGTGGATTCTCTCTGAAGCAGTTAACCTGG  
GGTAAAAAGTGGCAACCAAGTATATGTTGCACCGAGTGATCGAGGGCACTAAAACGGATGGCTCTCAACC  
TTTTACTCCTTCAAGTATCTAATCAAGGAGGGCAACTGCTCTGCTCAGAGTGGCCTCGCATGGCAGGACT  
GTGACTTCAAGGACGCTGAGGAAGCCGCACTGGAGAATGCACAGCAACTGTGGGAAAAGAGAAAAATGA  
ATTCTTCATAGTACCCAGACCTGCAAGATTGCTCCAAGTAAGGCCCCCACTACTGAAAGCCTATTTCCCC  
TGTATTGGTTGTGTGCATGCCATATCGACAGATAGTCCAGACCTGGAGCCTGTTCTGAAACTCCATCG  
AACATTTCAACAACAACACAGATCACAGCCACCTTTACTCTCAGAAAAGTAAAAAGTGCCACAGACA  
GGTGGTGGCTGGCCTGAATTTTGACATTACCTACACAATTGTGCAACAAATTGTTCAAAGGAGCGTTTT  
CCTTCCCTCCATGGAGACTGCGTGGCCCTTCCCAATGGTATGATGGTGAATGTAGAGGAAATCTCTTCA  
TGGATATTAATAACAAAATTGCCAACTTCTCACAGAGCTGTACCCTTTATTACAGGAGATGATTTGGTAGA  
AGCGCTTCCCAAGCCTTGCCCTGGCTGCCCCAGGGACATACCTGTAGACAGCCCAGAGCTGAAGGAGGTG  
CTTGGTCAATCCATTGCACAGCTAAATGCAGAGAATGACCATCCTTTCTATTACAAGATTGACACCGTGA  
AAAAAGCAACATCACAGGTGGTAGCAGGAACTAAATATGTTATTGAGTTCATAGCCAGAGAAACCAATG  
CTCCAAGGAAAGTAAACACAGAGCTGGCAGAAGATTGTGAGATCAAGCACCTTGGACAAAGTCTCGACTGC  
AATGCTAACGTGTACATGAGACCTTGGGAGAACAAGTCGTCCCGACTGTGAAATGCCAAGCATTAGATA  
TGACTGAAATGGCAAGAAGGCCCTCAGGTTTTTCTCCTTCCGGAGTGTACAGTACAAGAAAACAAAAGA  
AGGAAGAAGTGTAAAGTCCACCCTACATTGCCAGGGAGCAAGAAGAGAGGGATGCAGAAAAGTGAACAAGGA  
CCACACATGGGCATGGCTGGTTGCATGAAAAGCAAAATCAAGGCTAATAAGAATCACCGTGGTCATAAGC  
ATGGGCATGACCATGGCCATTGGTCCCCAAGGAGGCATGGCCTTGGTCATGGACACCAGAAACCACATGG  
CCTTGGTCATGGACATCAACTTAACTTGATTATCTTAGACACCAAGGGAAGATGGTATGACCATACA  
CACACAGTGGGACATGGTCATGGTCACGGACATGGTCATGGTCATGGACATGGACATGGTCATGGTCACG  
GTCATGGTCATGGTCATGGACATGGTCATGGTAAGCATACAAATAAAGACAAAAACAGTGTAAAGCAGAC  
TACCCAGAGAACAGAGTCTTTGGCAAGCTCTTCTGAATACAGTACTACATCTACACAGATGCAGGGGAGG  
ACAGAAGGGCCCACCTTGACCCCTCCCCGAGCCAGCCAAGTGTACCTCTTCTGGTTTTCAGGACTCAG  
ATTCATTGAAGATGTGGTAGCTACCACCCACCATATGACACTGGGGCCATGATGATTTGATCCCTGA  
TATCCATGTACAGCCAGATAGCCTTTCATTTAAGCTGATATCTGATTTTCCAGAAGCAACTTCCCCAAA  
TGTCTGGGCGCCATGGAAGCCAGCTAGCTGGAAGGATCCAACACAGAAAACAAGAATTTTCTGATT  
TTGATCTCCTTGATGCTCTTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



## Plasmid Map:



ACCN: NM\_001102411

ORF Size: 1983 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001102411.1](#), [NP\\_001095881.1](#)

RefSeq Size: 2299 bp

RefSeq ORF: 1986 bp

Locus ID: 16644

UniProt ID: [O08677](#)

Cytogenetics: 16 B1

MW: 73.6 kDa

**Gene Summary:** (1) Kininogens are inhibitors of thiol proteases; (2) HMW-kininogen plays an important role in blood coagulation by helping to position optimally prekallikrein and factor XI next to factor XII; (3) HMW-kininogen inhibits the thrombin- and plasmin-induced aggregation of thrombocytes; (4) the active peptide bradykinin that is released from HMW-kininogen shows a variety of physiological effects: (4A) influence in smooth muscle contraction, (4B) induction of hypotension, (4C) natriuresis and diuresis, (4D) decrease in blood glucose level, (4E) it is a mediator of inflammation and causes (4E1) increase in vascular permeability, (4E2) stimulation of nociceptors (4E3) release of other mediators of inflammation (e.g. prostaglandins), (4F) it has a cardioprotective effect (directly via bradykinin action, indirectly via endothelium-derived relaxing factor action); (5) LMW-kininogen inhibits the aggregation of thrombocytes; (6) LMW-kininogen is in contrast to HMW-kininogen not involved in blood clotting (By similarity).[UniProtKB/Swiss-Prot Function]