

Product datasheet for **MR217208**

Galnt6 (NM_001161767) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Galnt6 (NM_001161767) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Galnt6
Synonyms:	4632410F13; AW047994; GalNAc-T6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR217208 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGAGGTTGCTACGAAGCGCCACATGTCCTGCGCCTCGCCATGCTGGGACGCTCTTCATGCTCTTCC
 TCTTCATTCGGCAGAAGGATGTAAGCAACCAGGAGCAGGCCATGGAGAAACCGTGGCTGAAGTCCCTGGC
 TGGCCAGAAGGACCAGGTCCTGGACTTCATGCTAGGTGCTGTGAACAACATTAGGACGTCATGCCTAAG
 CTGCAGATCCGGGCTCCAGAACCGCCTCAGACTCTGGTCTCCACAAACCACTCTGCCTGCCCGGTTCT
 ACACCCCTGCTGAGCTGAAACCTTCTGGGACCGCCACCTCAAGATCCCAACAGCCCTGGGGCAGATGG
 AAAAGCATTTCAGAAGAAAGAATGGACCAATCTAGAGACCAAAGAAAAGGAAGGCTATAAGAAACAC
 TGCTTTAATGCATTTGCCAGTGACCGCATCTCCTTGACAGAGTCCCTGGGGCCAGACACCCGGCCTCCTG
 AGTGTGTTGACCAGAAGTCCGGCGCTGCCGCCACTGCCACCACCAGCGTGATCATTGTGTTCCACAA
 CGAAGCCTGGTCCACTCTGCTGCGCACAGTATACAGCGTCTGCACACCAGCCCTGCCATCTTGCTGAAG
 GAAATCATCCTCGTGGATGACGCCAGCAGATGAGCACCTGAAGGAGAGGCTGGAGCAGTACGTTACAGC
 AGCTGCAGATCGTGCGCGTGGTTCGGCAGCGGGAGCGGAAGGGGCTCATCACAGCCAGGCTGCTGGGGG
 CAGCGTGGCGCAGGCAGAGGTGCTCACGTTCTGGACGCCACTGTGAGTGTTCACGGCTGGCTGGAG
 CCCCTCCTGGCTCGAATTGCTGAAGACAAGACAGCGGTGGTGAGCCAGACATCGTCACCATCGACCTTA
 ATACCTTCCAATTCTCCAGACAGTACAGAGAGGCAAAGCCACAGCCGGGGCAACTTCGACTGGAGCCT
 GACCTTCGGCTGGGAAATGCTTCCTGAGCAGGAGAAGCAGAGACGCAAGGATGAGACCTACCCATCAAG
 TCCCGACGTTTCCGGCGGCCCTCTTCTCCATCTCCAAGGCCACTTTGAGCACATCGGTACCTATGATA
 ACCAGATGGAGATCTGGGGAGGGGAGAATGTGAAATGTCTTCCGGGTGTGGCAGTGTGGGGCCAGCT
 GGAAATAATCCCCTGCTCGGTGGTGGCCACGTGTTCCGTACCAAGAGTCCCACACCTTCCCAAGGGC
 ACCAGTGTATTGCGCGCAATCAGGTGCGCCTGGCTGAGGTCTGGATGGACGACTATAAGAGATTTTCT
 ACAGGAGAAATCTGCAGGCAGCAAAGATGGTCCAGGAGAACTTCGGTGACATCTCGGAGCGACTTCG
 GCTCAGGGAGCAGCTACGATGCCACAATCTCCTGGTACCTGCACAATGTCTACCCAGAGATGTTTGT
 CCCGACCTGAACCCACCTTCTATGGAGCCATCAAGAATCTTGGACCAATCAGTGTGGAGCTGGGCG
 AGAACAACCGGGAGGAAAGCCTCTCATCATGTATGTCTGCCACAACCTCGGCGCAATCAGTACTTTGA
 GTACACATCACAGAGGACCTCCGCCACAACATGGAAAGCAGCTGTGTCTACATGCCAGTGGGAGCAG
 CTGGCCTTAGGAGCTGTGAGTTCGTTGGCAAAAACAGCAGAGTGCCCAAGGACGAGGAGTGGGAATTGA
 CCCAGGATCAACTCATCAGAAACTCAGGATCTGGTACCTGCCTGACATCCCAGGACAAAAGCCAGCCAT
 GGCCCCGTGAATCCCAGGGATCCCTACCAGCTCTGGCTGTTTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR217208 protein sequence
 Red=Cloning site Green=Tags(s)

MRLLRRRHMSLR LAMLGSVFMFLFIRQKDVSNQEQAMEKPWLKSLAGQKQVLD FMLGAVNNIRDVMPK
 LQIRAPEPPQTLVSTNHSCLPGFYTPAELKPFWDRPPQDPNSPGADGKAFQKKEWTNLETKEEGYKHH
 CFNAFASDRISLQSLGPDTRPPECVDQKFRRCPLPTTSVIVFHNEAWSTLLRVTYVSLHTSPAILLK
 EIIILVDDASTDEHLKERLEQYVQQLQIVRVVRQERKGLITARLLGASVAQAEVLTFLDAHCECFHWLE
 PLLARIAEDKTA VVSPDIVTIDLNTFQFSRPVQRGKAHSRGNFDWSLTFGWEMLPEHEKQRRKDETYPIK
 SPTFAGGLFSISKAYFEHIGTYDNQMEIWGGENVEMSFVWQCGGLEIIPCSVVGHVFRKSPHTFPKG
 TSVIARNQVRLAEVWMDYKIFYRRNLQAAKMVQENNFGDISERLRLREQLRCHNFSWYLNHVYPEMFV
 PDLNPTFYGAIKNLGTNQCLDVGENNRGGKPLIMYVCHNLGGNQYFEYTSQRDLRHNIQKQLCLHASGST
 LGLRSCQFVGKNSRVPKDEEWELTQDQLIRNSGSGTCLTSQDKKPAMAPCNRPDPYQLWLFV

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001161767

ORF Size: 1869 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001161767.1](#), [NP_001155239.1](#)

RefSeq Size: 5715 bp

RefSeq ORF: 1869 bp

Locus ID: 207839


UniProt ID: [Q8C7U7](#)

Cytogenetics: 15 F1

MW: 71.5 kDa

Gene Summary: Catalyzes the initial reaction in O-linked oligosaccharide biosynthesis, the transfer of an N-acetyl-D-galactosamine residue to a serine or threonine residue on the protein receptor. May participate in synthesis of oncofetal fibronectin. Has activity toward Muc1a, Muc2, EA2 and fibronectin peptides (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR217208