

Product datasheet for **MR217263L4V**

Arid3c (NM_001017362) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arid3c (NM_001017362) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arid3c
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001017362
ORF Size:	1230 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217263).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001017362.2 , NP_001017362.1
RefSeq Size:	1751 bp
RefSeq ORF:	1230 bp
Locus ID:	550619
UniProt ID:	A6PWW5
Cytogenetics:	4 A5
Gene Summary:	This gene is a member of the ARID (AT-rich interaction domain) family of proteins. The ARID domain is a helix-turn-helix motif-based DNA-binding domain. ARID family members have roles in embryonic patterning, cell lineage gene regulation, cell cycle control, transcriptional regulation and possibly in chromatin structure modification. [provided by RefSeq, Jul 2008]



[View online »](#)