

Product datasheet for **MR217277L3V**

Adam33 (NM_033615) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adam33 (NM_033615) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adam33
Synonyms:	Adaml
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033615
ORF Size:	2391 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217277).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033615.2 , NP_291093.2
RefSeq Size:	3165 bp
RefSeq ORF:	2394 bp
Locus ID:	110751
UniProt ID:	Q923W9
Cytogenetics:	2 63.26 cM



[View online »](#)

Gene Summary:

This gene encodes a member of a disintegrin and metalloprotease (ADAM) family of endoproteases that play important roles in various biological processes including cell signaling, adhesion and migration. This gene is widely expressed, most highly in the adult brain, heart, kidney, lung and testis. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional metalloprotease enzyme. Alternative splicing results in multiple transcript variants encoding different isoforms, some of which may undergo similar processing. [provided by RefSeq, May 2016]