

Product datasheet for **MR217341L3V**

Dctd (NM_001161516) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dctd (NM_001161516) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dctd
Synonyms:	6030466N05Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001161516
ORF Size:	552 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001161516.1 , NP_001154988.1
RefSeq Size:	1877 bp
RefSeq ORF:	555 bp
Locus ID:	320685
UniProt ID:	Q8K2D6
Cytogenetics:	8 B1.1
Gene Summary:	Supplies the nucleotide substrate for thymidylate synthetase.[UniProtKB/Swiss-Prot Function]



[View online »](#)