

Product datasheet for MR217348L3

Atp13a2 (NM_001164366) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Atp13a2 (NM_001164366) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Atp13a2
Synonyms:	1110012E06Rik; AA589443
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217348).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001164366
ORF Size:	3345 bp



[View online >](#)

