

Product datasheet for **MR217426L3V**

Abca6 (NM_147218) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Abca6 (NM_147218) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abca6
Synonyms:	6330565N06Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_147218
ORF Size:	4872 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217426).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_147218.2 , NP_671751.2
RefSeq Size:	5240 bp
RefSeq ORF:	4875 bp
Locus ID:	76184
UniProt ID:	Q8K441
Cytogenetics:	11 E1
Gene Summary:	Probable transporter which may play a role in macrophage lipid homeostasis. [UniProtKB/Swiss-Prot Function]



[View online »](#)