

Product datasheet for **MR217470L4V**

Ube2z (NM_172300) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2z (NM_172300) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2z
Synonyms:	AW049199; C030047H17Rik; D11Moh35
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172300
ORF Size:	1068 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217470).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172300.3 , NP_758504.3
RefSeq Size:	3974 bp
RefSeq ORF:	1071 bp
Locus ID:	268470
UniProt ID:	Q3UE37
Cytogenetics:	11 59.28 cM
Gene Summary:	Catalyzes the covalent attachment of ubiquitin to other proteins. Specific substrate for UBA6, not charged with ubiquitin by UBE1. May be involved in apoptosis regulation. [UniProtKB/Swiss-Prot Function]



[View online »](#)