

Product datasheet for MR217535

Trappc5 (NM_025701) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Trappc5 (NM_025701) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Trappc5
Synonyms: 4021401A16Rik; TRS31
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR217535 representing NM_025701
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGCGGTTTCACACGCGGGAAGTCGGCGCTGCTAGAGCGCGCACTGGTGC GGCCCGCACCCGAGG
 TGAGCCTGAGCGCCTTTGCCCTTCTCTTCCGAGCTGGTGCAGCACTGTCAGAGCCGCTTTTCTCTGT
 GGCCGAGCTTCAAGCACGCTTGGCTGCCCTGGCCGACAGGTGGCGCCGAGTCTAGATGCCCTGGT
 GCCCGTGAGAAGGGTCCCGCGTGAGACTAAGGTGCTGGCGCATTGCTCTTCGTCAAGGGTGC GGGT
 GGAAGGCACTGTTCCGCAAGGAGGCGGACAAGCTGGAGCAGGCCAATGATGACGCCCGCACTTTTACAT
 CATCGAGAGGGAGCCGCTCATTAAACACCTACATCTCCGTGCCAAGGAGAACAGCACACTGAAGTGTGCC
 AGCTTACAGCGGGCATCGTGGAGGCTGTGCTCACGCACAGCGGCTTCCCTGCCAAGGTTACGGCGCACT
 GGCACAAGGGGACCACACTCATGATCAAATTTGAAGAGGCTGTGATTGCCCGAGATCGGGCCCTGGAGGG
 ACGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR217535 representing NM_025701
 Red=Cloning site Green=Tags(s)

MEARFTRGKSALLERLVRPRTEVSLSAFALLFSELVQHCQSRVFSVAELQARLALGRQVGARVLDALV
 AREKGARRETKVLGALLFVKGAVWKALFGKEADKLEQANDDARTFYIIEREPLINTYISVPKENSTLNCA
 SFTAGIVEAVLTHSGFPAKVTAHWHKGTTLMIKFEEAVIARDRALEGR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Chromatograms: https://cdn.origene.com/chromatograms/mm9059_e04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_025701

ORF Size: 564 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025701.4](#), [NP_079977.3](#)


RefSeq Size: 2072 bp

RefSeq ORF: 567 bp

Locus ID: 66682

UniProt ID: [Q9CQA1](#)
 Cytogenetics: 8 1.92 cM
 MW: 20.8 kDa
 Gene Summary: May play a role in vesicular transport from endoplasmic reticulum to Golgi.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR217535