

Product datasheet for **MR217544L3V**

Csnk1g3 (NM_152809) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Csnk1g3 (NM_152809) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Csnk1g3
Synonyms:	3300002K07Rik; C330049O21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_152809
ORF Size:	1272 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217544).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152809.2 , NP_690022.2
RefSeq Size:	4394 bp
RefSeq ORF:	1275 bp
Locus ID:	70425
UniProt ID:	Q8C4X2
Cytogenetics:	18 D1



[View online »](#)

Gene Summary:

Serine/threonine-protein kinase. Casein kinases are operationally defined by their preferential utilization of acidic proteins such as caseins as substrates. It can phosphorylate a large number of proteins. Participates in Wnt signaling (By similarity). Regulates fast synaptic transmission mediated by glutamate.[UniProtKB/Swiss-Prot Function]