

## Product datasheet for **MR217555L3V**

### Caprin1 (NM\_001111292) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Caprin1 (NM_001111292) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Caprin1
Synonyms:	AL022980; Caprin-1; Gpiap; Gpiap1; Mmgpip137; P137; rng105
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001111292
ORF Size:	2076 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217555).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001111292.1</a> , <a href="#">NP_001104762.1</a>
RefSeq Size:	3508 bp
RefSeq ORF:	2079 bp
Locus ID:	53872
UniProt ID:	<a href="#">Q60865</a>
Cytogenetics:	2 E2



[View online »](#)

**Gene Summary:**

May regulate the transport and translation of mRNAs of proteins involved in synaptic plasticity in neurons and cell proliferation and migration in multiple cell types (PubMed:20516077). Binds directly and selectively to MYC and CCND2 RNAs. In neuronal cells, directly binds to several mRNAs associated with RNA granules, including BDNF, CAMK2A, CREB1, MAP2, NTRK2 mRNAs, as well as to GRIN1 and KPNB1 mRNAs, but not to rRNAs (By similarity).[UniProtKB/Swiss-Prot Function]