

Product datasheet for **MR217562**

Nip7 (NM_001164472) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nip7 (NM_001164472) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nip7
Synonyms: 1110017C15Rik; 6330509M23Rik; AA408773; AA410017; PEachy
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR217562 representing NM_001164472
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGGCCCTTGACTGAAGAAGAGACCCGGTCATGTTTGAGAAGATCGGAAGTACATCGGAGAGAACC
TCAACTCCTAGTGGACAGACCCGACGGCACGTACTGTTCCGGCTGCACAACGACCGAGTGTACTATGT
GAGTGAGATGATGCTGAAGCTGGCCGCAACATCTCCGGAGACAAGCTCGTGTGCTGGGACGTGTTTT
GGGAAGTTTACCAAGACCCACAAGTTCGCTTGCACGTGACAGCTCTGGATTACCTGGCCCCCTATGCCA
AGGGTTTTGGGGTGGCAGCAAAGTCTACCCAAGACTGCCGAAAAGTGGACCCCATGGCAATTGTGGTGT
CCACCAAGCAGACATTGGGGAGTACGTACGCCATGAAGAGACTGACT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR217562 representing NM_001164472
Red=Cloning site **Green**=Tags(s)

MRPLTEEETRMFEKIAKYIGENLQLLVDRPDGTYCFRLHNDRVYVSEMLKLAANISGDKLVSLGTCF
GKFTKTHKFRHLHTALDYLAAPYAKGFGVAAKSTQDCRKVDPMAIVVFHQADIGEYVRHEETLT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001164472

ORF Size:

399 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001164472.1</u> , <u>NP_001157944.1</u>
RefSeq Size:	2895 bp
RefSeq ORF:	402 bp
Locus ID:	66164
Cytogenetics:	8 D3
MW:	15.7 kDa
Gene Summary:	Required for proper 34S pre-rRNA processing and 60S ribosome subunit assembly. [UniProtKB/Swiss-Prot Function]