

Product datasheet for MR217584

Chn1 (NM_001166603) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Chn1 (NM_001166603) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Chn1
Synonyms: 0610007I19Rik; 0710001E19Rik; 1700112L09Rik; 2900046J01Rik; AI413815; ARHGAP2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR217584 representing NM_001166603
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTCGACCCTGTTTCGATACAGATGAATATAGACCTCCTGTTTGGAAATCTTACTTATACCAGCTGC
 AGCAGGAAGCCCCTCACCTCGAAGGGTCACCTGCACTTGTGAGGTAGAAAACAGACCAAAGTATTATGG
 AAGAGAGTATCATGGCATGATCTCTAGAGAAGAGACTGACCAGCTCCTGAGTGTGGCTGAGGGGAGCTAC
 CTCATCCGTGAGAGCCAGCGCCAGCCAGGAACGTACACTTTGGCTTTAAGATTTGGAAGTCAGACCAGAA
 ACTTCAGGCTGTACTACGATGGAAGCACTTTGTTGGGAGAAAACGCTTTGAGTCCATCCACGATCTGGT
 GACTGATGGCTTGATTACTCTCTATATTGAAACCAAGGCAGCAGAATACATTGCCAAGATGACGATAAAC
 CCAATTTATGAGCACATAGGATACACAACCTTAAACAGAGAGCCAGCATACAAAACAGCACATGGCAGTCC
 TGAAAGAGACACATGATGAGAAAAGAGGCTACAGGCCAGGATGGGGTATCAGAGAAAAGGTTGACATCACT
 TGTTAGAAGAGCAACTCTGAAAGAAAATGAACAAATTCAAAAATGAAAAGGTTCACAATTTCAAGGTG
 CATACGTTCCGAGGGCCACACTGGTGTGAATACTGTGCCAACTTCATGTGGGGCCTCATTGCTCAGGGAG
 TGAAATGTGCAGGTCTTAATTCTGAAGGACTCTACCGAGTGTGAGGATTTAGTGACCTGATTGAAGATGT
 CAAGATGGCTTTTGATAGAGATGGTGAGAAGCCGATATTTCTGTGAACATGTATGAGGACATCAACATT
 ATCACTGGTGCCTTAACTGTACTTCAGGGATCTGCCAATTCCTCATCACATACGATGCCTACCCCA
 AGTTTCATTGAGTCTGCCAAAATATGGACCCTGACGAGCAATTGGAGACCCTTCACGAAGCACTGAGATC
 GCTGCCGCTGCCACTGCGAGACGCTCCGGTACCTCATGGCGCATCTCAAGAGAGTGACCCCTCATGAG
 AAGGAGAATCTGATGAGTGCAGAGAACCTTGGGATCGTGTGGACCAACCCTCATGAGATCCCCAGAGC
 TCGACCCCATGGCCGCCCTGAACGACATACGCTATCAGAGACTGGTGGTGGAGCTGCTTATCAAAAACGA
 AGACATTTTATT

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >MR217584 representing NM_001166603
Red=Cloning site Green=Tags(s)

MALTLFDTDEYRPPVWKSYLEQLQEQEAPHPRVTCTCEVENRPKYYGREYHGMISREETDQLLSVAEGSY
 LIRESQRQPGTYTLALRFGSQRNFRLLYYDGKHFVGEKRFESIHDLVTDGLITLYIETKAAEYIAKMTIN
 PIYEHIGYTTLNREPAYKQHMVAVLKETHDEKEATGQDGVSEKRLTSLVRRATLKENEQIPKYEKVHNFKV
 HTFRGPHWCEYCANFMWGLIAQGVKCAGLNSEGLYRVSGFSDLIEDVKMAFDRDGEKADISVNMVEDINI
 ITGALKLYFRDLPIPLITYDAYPKFIESAKIMDPDEQLETLHEALRSLPPAHCETLRYLMAHLKRVTLHE
 KENLMSAENLGI VFGPTLMRSPELDPMALNDIRYQRLVVELLIKNEDILF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001166603

ORF Size: 1203 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001166603.1, NP_001160075.1

RefSeq Size: 3876 bp

RefSeq ORF: 1206 bp


Locus ID: 108699

Cytogenetics: 2 C3

MW: 47 kDa

Gene Summary: GTPase-activating protein for p21-rac and a phorbol ester receptor. May play an important role in neuronal signal-transduction mechanisms (By similarity). Involved in the assembly of neuronal locomotor circuits as a direct effector of EPHA4 in axon guidance.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR217584