

Product datasheet for **MR217585L4V**

Tead3 (NM_001098226) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tead3 (NM_001098226) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tead3
Synonyms:	DTEF-1; ETRF-; ETRF-1; Tcf13r; Tcf13r2; TEAD-; TEAD-3; TEF-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001098226
ORF Size:	1395 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217585).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001098226.3 , NP_001091696.2
RefSeq Size:	2448 bp
RefSeq ORF:	1398 bp
Locus ID:	21678
UniProt ID:	P70210
Cytogenetics:	17 14.66 cM



[View online »](#)

Gene Summary:

This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and thought to play a role in placental gene regulation and development. Alternative splicing, and alternate use of an upstream AUG translation initiation codon, and an in-frame downstream non-AUG (AUA) codon, results in 2 isoforms. [provided by RefSeq, Jul 2008]