

Product datasheet for **MR217593L4V**

AtI3 (NM_001163505) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AtI3 (NM_001163505) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	AtI3
Synonyms:	4633402C03Rik; 5730596K20Rik; AI465397; AW228836
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001163505
ORF Size:	1608 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217593).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001163505.1 , NP_001156977.1
RefSeq Size:	6485 bp
RefSeq ORF:	1611 bp
Locus ID:	109168
UniProt ID:	Q91YH5
Cytogenetics:	19 A
Gene Summary:	GTPase tethering membranes through formation of trans-homooligomers and mediating homotypic fusion of endoplasmic reticulum membranes. Functions in endoplasmic reticulum tubular network biogenesis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)