

## Product datasheet for **MR217597L4V**

### Hecw2 (NM\_001001883) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hecw2 (NM_001001883) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hecw2
Synonyms:	A730039N16Rik; D030049F17Rik; Gm971; mKIAA1301; Nedl2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001001883
ORF Size:	4734 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217597).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001001883.3</a> , <a href="#">NP_001001883.1</a>
RefSeq Size:	11179 bp
RefSeq ORF:	4737 bp
Locus ID:	329152
UniProt ID:	<a href="#">Q6I6G8</a>
Cytogenetics:	1 C1.1
Gene Summary:	E3 ubiquitin-protein ligase that mediates ubiquitination of TP73. Acts to stabilize TP73 and enhance activation of transcription by TP73. Involved in the regulation of mitotic metaphase/anaphase transition.[UniProtKB/Swiss-Prot Function]



[View online »](#)