

## Product datasheet for **MR217626**

### **Kbtbd11 (NM\_029116) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Kbtbd11 (NM_029116) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Kbtbd11
Synonyms:	2900016B01Rik; 4930465M17Rik; mKIAA0711
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR217626 representing NM\_029116  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGAGAACTCGGTGGCGCTTTCGTCCTACTCGGGACTGAGCCTAGGACTCCTGGAGAGGACAGCC  
TCCTCTGCCAGCCGAGGAAGAGGGCGCTGCGTCCACAGCGCAGACACCCTGCAGTCTCAGCGCCTCATT  
GTGCTTACAGTCCGGGGACGACTCCCGCCGAGTCTCGCGCTCTGCCGCGGAGGGCTCCGAAGCCTCT  
CCGCTTCTCTCCGACGACCTGCGTGTGGTGGAGACGAGTGGGATGTCAGCAGCGCCGCTCCCGAG  
AGTCCCGGAGGAGTGCAGCGCCCGGAAGAGCCCGCTCGCGGAGGATCCCCCTCCGACACGAGCA  
CGCGCTCCCGTGAACCTGAGTCTTTGGACGAGCTTGGGAGCCCGTGCCTGTGCCGCCGGGGTGGG  
TCGGTGCACGGAGAACCAGCTGGTATCGAGGTGGCCGGTCCCGACTACGGGCGCATAAGGCGGTGC  
TGGCCGACGCGAGGACTACTCCGCGCGCGCATCCCGGACGTGCTGCGGGTCAAGCGTGAGCTT  
CACTGCCCTGCGCTACTGCTGGCCGACGACATACAGCGGGCGCATGGCGGGCTGCGGCCGACAACGTG  
GCCGAGGTGGTGGCCGGCGCGCCGCCTCCAGCTGCCCGGCGCTGCGCAGCGCGCCACCGAGGCCATGG  
CGCCGACGCTGAGCCTGGACAACCTGCTACGAGGTGCTGAGCGCAGGCAAGCGGCAGCGGGTACCCGAGCT  
GCGCGACGCCGCTACCGCTTATGAGCGACCACTACCTGGAGGTGCTGCGCGAGCCCGCCTATTCGGC  
CGCCTGTCCGGTGCAGAACGCGATCTGCTGCTGCGGCGCCGCTGTGTACCGGCCGCGCTTGTCTGCTGG  
CCGCTGCGCTCGGGACAACCGGGGAGCGCTTGGCAGTCGACCGCAGAGCCCGTCCGGGGACGAGATC  
GCGCGGGGACGCGGCCGTTTATTGCTACCAGGCGGAGCCGGGGAGTGGCGGAGCTGACCCGGCTGCC  
GAGGGCGCGCCAGCGCGGGCTGCGGCCCTGCGTGTCTTCAACTACCTTCTCTGGCGGGCGGTGG  
CTCCTGCGGGTCCCGATGGCCGCGCGCCCGTCCGACCAGGTGATTGCTACAACCCGGTACCCGACTC  
CTGGAGTACTGTGCGCCGCTGCCAGGACGCTCGCAGGTGCAGCTGCTGGCCCTGGACGGCCACCTG  
TATGCCGTGGGTGGCAGTGCCTGCTCAGCGTGGAGCGCTACGACCCGCGAGCTGACCGATGGACGGCGG  
TGGCCCGCTTCCCGTGGCGCCTTCGAGTGGCCACGAGGCTGCCACCTGTAACGGAGAGATTTATGT  
GTCCGGGGGCTCACTTCTACCGCTGCTCAAGTACGACCCCGCGCGACGAGTGGCAGGAGTGCCTCA  
TGCAGCAGCAGCCGAGAGCGCTTGCAGCATGGTGGCCCTCGACGGCTTCTCTATCGCTTTGACCTGT  
GTGGGAGCCGTGGGAAGCTCAGGCTGCAGTCGGTTCAGGTGGCGGGTACGCTTTTCCGCTACCATTG  
CCTGGCCAAGCAATGGAGCAATGTGCTGTACACCTGAGGCCCTGGCGACCCGAGGCTCCAGCCC  
TTTCGTTGTGTTGCCCTGGATGGTACCATCTACTGTGTGAGCCGCGCGGACCTGGCGTTTCGTGCCTT  
CTCAGGACACCGAGGCTGGCAGCGATATGGACCTGGTGGCAGCTTTGAGCCCGAGCCCCTGGGATCCCC  
CCTGGATGTCCGGGGCTACTCTTCCATTTGACTCAACCTGCCTGAGAAGCCGACCGCGGGGAGCAG  
GGTGCAGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTAA

**Protein Sequence:** >MR217626 representing NM\_029116  
 Red=Cloning site Green=Tags(s)


MENSVAPFVLYSGTEPRTPGEDSLPLPAEEEGAASTAQTPCSLSASLCFSSGDDSPQSRASAAEGSEAS  
 PPSLRSDLRVVETQWDVSSAASPESPEECARPEEPASPEDPPSRHEHARPVELESDELGEPVPPPGVG  
 SVHGEPDLVIEVAGRRLRAHKAVLAARSDYFRARASRDVLRVQGVSF TALRLLLADAYSGRMAGVRPDNV  
 AEVVAGARRLLQPGAAQRATEAMAPQLSLDNCYEVLSAGKRQRL TELRDAAYRFMSDHYLEVLRPAVFG  
 RLSGAERDLLLRRLCTGRACLLAAALGTTGERSGRPQSPSGDAESRGDAAVYCYQAEAGEWRELTRLP  
 EGAPARGCGLCVLFNYLFLAGGVAPAGPDGRARPSDVVYCNPVTDWSWTVRPLRQARSQVQLLALDGH  
 YAVGGECLLSVERYDPRADRWTAVAPLPRGAFAVAHEAATCNGEIYVSGGSLFYRLLLKYDPRRDEWQEQCP  
 CSSSRERSADMVALDGLYRFDLCGSRGEAQAAVGSGGGVSVFRYHCLAKQWSQCAVHLRPPGAPAGLQP  
 FRCVALDGTIYCVSRAGTWRFPVPSQDTEAGSDMGPGGSFEPPEPLGSPLDVRGVLFPFVLNLPEKPDGRGEQ  
 GAV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_029116

**ORF Size:** 1899 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_029116.2](#), [NP\\_083392.1](#)

**RefSeq Size:** 6996 bp

**RefSeq ORF:** 1902 bp


**Locus ID:** 74901

**UniProt ID:** [Q8BNW9](#)

**Cytogenetics:** 8 A1.1

**MW:** 67.9 kDa

Product images:



Circular map for MR217626