

Product datasheet for **MR217633L3V**

Tmem235 (NM_001085535) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem235 (NM_001085535) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem235
Synonyms:	Cldn27; Gm12581; OTTMUSG00000007589
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001085535
ORF Size:	633 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217633).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001085535.1 , NP_001079004.1
RefSeq Size:	1820 bp
RefSeq ORF:	636 bp
Locus ID:	546519
UniProt ID:	B1AQL3
Cytogenetics:	11 E2



[View online »](#)