

## Product datasheet for **MR217667L4V**

### Asb8 (NM\_030121) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Asb8 (NM_030121) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Asb8
Synonyms:	4930539L19Rik; AI788835; C430011H06Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030121
ORF Size:	864 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217667).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030121.4</a> , <a href="#">NP_084397.3</a>
RefSeq Size:	2166 bp
RefSeq ORF:	867 bp
Locus ID:	78541
UniProt ID:	<a href="#">Q91ZT9</a>
Cytogenetics:	15 F1
Gene Summary:	May be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]



[View online »](#)