

Product datasheet for **MR217670L3V**

Fyco1 (NM_001110253) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fyco1 (NM_001110253) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fyco1
Synonyms:	2810409M01Rik; Mem2; RUFY3; ZFYVE7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001110253
ORF Size:	4314 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217670).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001110253.1 , NP_001103723.1
RefSeq Size:	7695 bp
RefSeq ORF:	4314 bp
Locus ID:	17281
UniProt ID:	Q8VDC1
Cytogenetics:	9 74.52 cM
Gene Summary:	May mediate microtubule plus end-directed vesicle transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)