

## Product datasheet for **MR217690L3V**

### Papss1 (NM\_011863) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Papss1 (NM_011863) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Papss1
Synonyms:	AI325286; Asapk; SK1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011863
ORF Size:	1872 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217690).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011863.2</a> , <a href="#">NP_035993.1</a>
RefSeq Size:	2617 bp
RefSeq ORF:	1875 bp
Locus ID:	23971
UniProt ID:	<a href="#">Q60967</a>
Cytogenetics:	3 61.05 cM



[View online »](#)

**Gene Summary:**

Bifunctional enzyme with both ATP sulfurylase and APS kinase activity, which mediates two steps in the sulfate activation pathway. The first step is the transfer of a sulfate group to ATP to yield adenosine 5'-phosphosulfate (APS), and the second step is the transfer of a phosphate group from ATP to APS yielding 3'-phosphoadenylylsulfate (PAPS: activated sulfate donor used by sulfotransferase). In mammals, PAPS is the sole source of sulfate; APS appears to be only an intermediate in the sulfate-activation pathway (PubMed:7493984). Required for normal biosynthesis of sulfated L-selectin ligands in endothelial cells (By similarity). [UniProtKB/Swiss-Prot Function]