

Product datasheet for **MR217700**

Sestd1 (NM_175465) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sestd1 (NM_175465) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sestd1
Synonyms:	1500031J16Rik; A930010O20; Solo/Trio8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR217700 representing NM_175465
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGGAGGCTCGGTGATATTGCCATTCTGAAGAAGAAGCTGGCCTTCCTTTTCAGGAGAAAGGACAGAC
GCAGTGGACTCATCTTGACAATTCGCTGTGCCTTGAGCAGACAAGTATGGATGAGCTGAGTGCACTTT
GGATTATCTACTCAGCATTCCAAGTGAGAAATGAAGGCCAGAGGATTCCTGTGATCGTGGATGGCAGG
AAGTCTCAGTGAATGTGGTGAAGACAGTTGTCTTAATGCTGCAGAATGTTGTTCCAGCCGAGGTATCCC
TTGTTTGTGAGTGAAGCCAGATGAATTCTGGGATAAGAAAGTAACACACTTTTGTGTTTGGAAAGAAA
GGACAGACTTGGCTTTGAGGTTATTTAGTGTCCGCAATAAGTTGACGCGTTACATAGAGCCATGCCAG
TTGACAGAAGACTTTGGTGGGAGTCTCACCTATGACCACATGGACTGGCTGAATAAGAGGCTGGTCTTTG
AGAAGTTTACAAAGGAGTCTACATCATTGCTAGATGAAGTCTGCTGATTAACAATGGAAGTGATAAGGG
GAATGAGCAGGAGAAAGAAAGGCTGTGGATTTAACTTTCTCCATCCGTTGATCCTGAGACGGTCTCT
CAGACAGGGCATGAGTTGTTATCTGAATTACAGCAGCGTCGATTTAATGGCTCCGACGGAGGGTCTCGT
GGTCTCCTATGGATGACGAGCTGCTGGCACAGCCGAGGTTATGAAGTTGCTGGACTCACTCCGGGAGCA
GTACACGCGCTACCAGGAAGTCTGCAGGCAGCGCAGTAAGCGCACACAGCTGGAGGAGATCCAGCAGAAG
GTGATGCAGGTTGTGAAGTGGCTTGAAGGGCTGGATCTGAGCAACTGAGAGCACAGTGGGGCATTGGTG
ACTCGATCAGGGCCTCCCAGGCCCTGCAACAGAAGCATGAGGAGATCGAGAGCCAACACAGCGAGTGGTT
CGCAGTACGTGGAGCTTAACACAGCAGATTGCAGCACTTTAAACGCTGGAGACGAGGAAGATCTTGTG
GAGTTGAAGTCTCTACAGCAACAGCTTAGCGACGTTTGCATCGACAGGCCAGTCAGCTAGAATTTAGGC
AAAACCTCCTGCAGGCAGCTTTGGAATTTATGGTGTGGCCCAAGATTTGCTCAGCAGTTGGACGGCTT
ATTAGGGATGTTGTGCGTAGATGTAGCACCAGCTGATGGAGCATCGATTACAGAACTTTAAACTGCTT
GAAGAGAAGCTGAAAAGTGTGATGTAGGATTGCAAGGCTTGCAGAAAAAGGTCAAGGTCTCCTGGATC
AGATCTCCAACAGGCTTCTGGCCTATGGCAAGGATGTAAGTATTGAAAAAAGAAAACGTGGACCA
CATACAAGGAGTGTGGAAGACATGCAACTTAGGAAACAAAGATGTGAGGACATGGTAGATGTGCGGAGG
CTGAAGATGCTTCAGATGGTGCAGTTGTTAAGTGTGAAGAGGATGCTTCACAGGCAGTAGAATGGCTGA
GTGAAGTCTGGATGCCCTGCTGAAGACCCATATCAGGCTGGGTGACGATGCTCAGGAAACAAAGTTTTT
ACTGAAAAACACAGAAAATTTGTCGATGTTGCCAGAGTACTTACGACTATGGCAGACAGCTGCTGCAG
GCCACAGTTGACTGTGCCAGTCTCTGCGCTGCACTTCCCGGTCCTCAGGGGACACACTTCTCGACTGA
ACAGAGTGTGGAAGCAGTTTACAGTCGCATCTGAAGAGAGAGTGCACAGGTTGGAGATGGCTATTGCCTT
TCACTCCAATGCCGAAAAGATTTTGAAGACTGCCAGAAAGAGCCTGAAGCAATGAATGATGAGGAACAG
TTTGAGGAAATGGAAGCAATTGAAAAATCACTTCTGGATAGACTCACTATCCCCGTGGTTTACCCCGACG
GAACTGAACAATTTTTGGAAGTCCAAGTGACATGGCTTCTACTGCAGAGCACATCAGGGACAGAATGAA
ACTGGTCAGTCTCAAAGGCAGCAGCTGAGACATCCCGAATTGGTGACCACGGAGAGC

AGCGGACCGACGCGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR217700 representing NM_175465
 Red=Cloning site Green=Tags(s)


MEASVILPILKKKLAFLSGGKDRRSLILTIPLCLEQTSMDLSVTLDYLLSIPSEKCKARGFTVIVDGR
 KSQWNVVKTVVLMQNVVPAEVS LVCVVKPDEFWDKVKVTHFCFWKEKDRLGFEVILVLSANKLTRYIEPCQ
 LTEDFGGSLTYDHMDLWLNKRLVFEKFTKESTSLDELALINNGSDKGNEQEKERSVDLNF LPSVDPETVL
 QTGHELLSELQRRFNGSDGGVSWSPMDELLAQPOVMKLLDSLREQYTRYQEVCQRQRSKRTQLEEIQQK
 VMQVWNWLEGPGEQLRAQWIGDSIRASQALQKHEEIESQHSEWFAYVELNQIAALLNAGDEEDLV
 ELKSLQQQLSDVCYRQASQLEFRQNLQAALFEHGAQDLSQQLDGLLGMLCVDVAPADGASIQQLKLL
 EEKLSVDVGLQGLREKGGQLLDQISNQASWAYGKDVTIENKENVDHIQGYMEDMQLRKQRCEDMVDVRR
 LKMLQMVQLFKCEEDASQAVEWSEL LDALLKTHIRLGDDAQETKVLLEKHKRFVDVAQSTYDYGRLLQ
 ATVVLQSLRCTSRSSGDTL PRLNRVWKQFTVASEERVHREMAIAFHSNAEKILQDCPEEPEAMNDEEQ
 FEEIEAIGKSLDLRLTIPVVYPDGTEQYFGSPSDMASTAEHIRDRMKLVSLKRQQLRHPELVTTES

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9098_e07.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_175465


ORF Size: 2088 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM_175465.6](#), [NP_780674.1](#)
- RefSeq Size:** 8951 bp
- RefSeq ORF:** 2091 bp
- Locus ID:** 228071
- UniProt ID:** [Q80UK0](#)
- Cytogenetics:** 2 C3
- MW:** 79.4 kDa
- Gene Summary:** May act as the primary docking protein directing membrane turnover and assembly of the transient receptor potential channels TRPC4 and TRPC5. Binds phospholipids such as phosphatidylinositol monophosphates, phosphatidylinositol diphosphates (PIP2s) and phosphatidic acid, but not less polar lipids including phosphatidylcholine, phosphatidylserine, and phosphatidylinositol. The binding to PIP2s is calcium dependent. Might be involved in the plasma membrane localization of CTNNB1 (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR217700