

Product datasheet for **MR217730A1V**

Mouse *Dusp15* (NM_001159376) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	Dusp15
Synonyms:	AI851682; LMW-DSP10; T-DSP10
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR217730 representing NM_001159376
Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GGCGCGCC**

ATGGGGAATGGCATGACCAAGGTACTTCTGGACTCTACCTTGAACTTCATTGATGCCAAGACCCGG
ATCAGCTGGGCCGGAATAAGATCACACATATCATCTCTATCCAGTAATACCCAGCCTCTGTTGCAGGA
TATCAGCTACCTTCGAATCTCAGTGTCTGATACTCCTGAGGTACCCATCAAAAAGCACTTCAAAGAATGC
GTCCACTTTATCCACTCTGCCGCTCAACGGGGGTAAGTGCCTTGTGCACTGCTTTGCAGGTATATCTC
GAAGTACCACCATAGTGATCGCCTATGTGATGACGGTGAAGTGGCTGGCAGGAAGTGTGGAAGC
CATCAAAGCCTCCCGGCCATCGCAACCCAAACCCGGGCTTAGGCAGCAGCTCGAGGAGTTTGGCTGG
GCAAACCTCCAGAAGCTTCGCCGCGAGCTGGAGGAGCGCTTCGGGGAGATCCCGTTCCGCGACGAGGAAG
ACTTGCGCGCGCTGCTGCCTCTCTGCAGGCGCTGTGCCAGGGCTCGGCACATCGGCAGCGTCGCCAC
CACAGCTTCTCGGCCGCTGAGGGGACCTGCAGCGCCTGGTGCCGCGATCGCCGCGGGACTCACACCAG
CCGCTGCCGCTGTGGCGCGCTCAAGCAGACTTCTTCTGCCTCCCGCGGTGCCTGTGCGCAAGGGCG
GCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR217730 representing NM_001159376
Red=Cloning site Green=Tags(s)

MNGMTKVLPGLYLGNFIDAKDPDQLGRNKITHIISIHESPOPLLQDITYLRISVSDTPEVPIKKHFKEC
VHFHSCRLNNGNCLVHCFAGISRSTTIVIAVMTVTGLWQEVLEAIKASRPANPNPGRFQQLEEFGW
ANSQKLRRQLEERFGEIPFRDEEDLRALLPLCRRCRQGSATSAASATTASSAAEGLQLRVPSPRDSHQ
PLPLLARVKQTFCLPRCLSRKGGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001159376
ORF Size:	705 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001159376.1 , NP_001152848.1
RefSeq Size:	3121 bp
RefSeq ORF:	708 bp
Locus ID:	252864
UniProt ID:	Q8R4V2
Cytogenetics:	2 H1
MW:	26.2 kDa