

## Product datasheet for MR217750

### Commd6 (NM\_001168592) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Commd6 (NM\_001168592) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Commd6  
**Synonyms:** 1110059J08Rik; 1700063H17Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR217750 representing NM\_001168592  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGGAGTCTGGTTTCAGGGAGCCCGTCTGGATGCCAAGTCGGAGCTTATAGATTTTCAGTGGAAAC  
TGGGCATGGCTGTGAGCTCTGACAGCTGCAGATCTCTCAAGTATCCTTACGTGGCAGTGATGCTGAAGGT  
GGCAGATCACTCTGGCCAAGTAAGCAGCAAGTCCATCGAGATGACAATCCACAATTTCAGAATTTCTAC  
AAACAGTTCAAAGAAATTGCCGCTGTAATTGAGACTGTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR217750 representing NM\_001168592  
Red=Cloning site Green=Tags(s)  
MEESGFREPVLDAKSELIDFQWKLGMVSSDSCRSLKYPYVAVMLKVADHSGQVSSKSIEMTIPQFQNFY  
KQFKEIAAVIETV

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001168592

**ORF Size:** 249 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


**RefSeq:** [NM\\_001168592.1](#), [NP\\_001162064.1](#)

**RefSeq Size:** 892 bp

**RefSeq ORF:** 252 bp

**Locus ID:** 66200  
**Cytogenetics:** 14 E2.3  
**MW:** 9.9 kDa  
**Gene Summary:** May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes. Inhibits TNF-induced NFKB1 activation.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR217750