

Product datasheet for **MR217779**

Col4a3bp (NM_023420) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Col4a3bp (NM_023420) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Col4a3bp
Synonyms:	2810404O15Rik; 9230101K08Rik; AU016711; CERT; GPBP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR217779 representing NM_023420
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTCGGATAACACAGAGCTGGAACCTCGTCGGGCTCGGAGGAGGATCCGGAGACGGAGTCCGGGCCGCTG
 TGGAGCGCTCGGGGTCCTCAGCAAGTGGACAACTATATTCATGGATGGCAGGATCGTTGGGTAGTTTT
 GAAAAATAACTTTGAGTTACTACAACTCTGAAGATGAAACAGAATATGGCTGTAGGGATCCATCTGT
 CTTAGCAAGGCTGTGATCACGCCTCACGATTTTGATGAATGCCGGTTTGATATCAGTGTAAATGATAGTG
 TTTGGTACCTTCGAGCTCAGGACCCGGAGCACAGACAGCAATGGGTAGACGCCATTGAACAGCACAAAGAC
 TGAATCGGGATATGGATCTGAGTCCAGCTTGCCTAGACATGGCTCAATGGTGTCACTGGTGTCTGGAGCG
 AGTGGCTATTCTGCTACGTCACCTCTTCTTTCAAGAAAGGCCACAGTTTACGTGAGAACTGGCTGAAA
 TGGAGACATTTCCGGACATCCTGTGCCGCGAGTTGATACTCTCCAGAAGTACTTTGATGTCTGTGCTGA
 CGCTGTCTCCAAGGATGAGCTTCAGAGGGATAAAGTCGTAGAAGATGATGAAGATGACTTCCCTACAAC
 CGTTCTGATGGAGACTTTTTGCACAATACCAATGGTAATAAAGAAAAATTATTTCCACATGTAACACCAA
 AAGGAATTAATGGCATAGACTTTAAAGGGGAAGCAATAACTTTTTAAAGCAACTACTGCTGGAATCCTTGC
 TACACTTTCTCATTGTATTGAATTAATGGTAAAACGGGAAGAGAGCTGGCAAAAAAGACACGATAGGGAA
 GTGGAAAAGAGGAGACGAGTGGAGGAAGCGTACAAGAAATGTGATGGAAGAAGCTTAAAGAAGAAACCCGTT
 TCGGAGGGCCGGATTATGAAGAAGGTCCAAACAGTCTGATTAATGAGGAAGAGTTCTTTGATGCTGTTGA
 AGCTGCTCTTGACAGACAAGATAAAATAGAGGAACAGTACAGAGTGAAGGTCAGGTTACACTGGCCC
 ACATCATTGCCATCTGGAGACACCTTTTCTTCTGTCCGGACGCATAGATTTGTACAAAAGCCCTATAGTC
 GCTCTTCCCTCCATGTCTCCATTGATCTAGTCACTGCTGACGATGTTACAGATTGATGCAAAATGGCAACTG
 TGAAGAAATGGTACAGAACCACATGAACATTCATTACAGGATGTAGGTGGTGTGCAAAATGGCAACTG
 GTTGTGGAAGAAGGAGAAATGAAGGTATACAGAAGAGAAGTGAAGAAAAATGGAATTGTTCTGGATCCTT
 TGAAAGCTACTCATGCAGTTAAAGGTGTTACAGGACATGAGGTCTGCAATTACTTTGGAATGTTGATGT
 TCGCAATGACTGGGAACTACTATAGAAACTTTTATGTGGTGGAAACATTAGCTGATAATGCAATCATC
 GTTTATCAAACGCACAAGAGAGTATGGCCCGCTTCTCAGAGAGACGACTGTATCTTTCTGCTATTCGAA
 AGATCCCAGCCTTGACTGAAAATGATCCTGAACTTGGATAGTTTGTAAATTTTCTGTGGATCATGATAG
 TGCTCCTCTGAACAATCGATGTGTCGTGCCAAAATCAATATTGCTATGATTTGTCAAATTTAGTAAGC
 CCACCAGAGGGAGACCAGGAGATAAGCAGAGACAACATCTGTGCAAGATCACGTATGTAGCTAATGTGA
 ACCCAGGAGGATGGGCGCCAGCCTCGGTCTTAAGAGCAGTGGCAAGCGAGAATACCTAAGTTTCTAAA
 ACGTTTTACTTCTTATGTCCAAGAAAAACTGCAGGAAAACCAATTTTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR217779 representing NM_023420
 Red=Cloning site Green=Tags(s)

MSDNQSWNSSGSEEDPETESGPPVERCGVLSKWTNYIHGWQDRWVVLKNNTLSYYKSEDETEYGCGRGSIC
 LSKAVITPHDFDECRFDISVNDVWYLRQDPPEHRQWVDAIEQHKTESGYGSESSLRRHGMVSLVSGA
 SGYSATSTSSFKKGHSLREKLAEMETFRDILCRQVDTLQKYFDVCADAVSKDELQRDKVVEDDEDDFPPT
 RSDGDFLHNTNGNKEKLFPHVTPKINGIDFKGEAITFKATTAGILATLSHCIELMVKREESWQKRHDRE
 VEKRRRVEEAYKNVMEELKKKPRFGPDYEEGPNSLINEEEFFDAVEAALDRQDKIEEQSQSEKVRHLHP
 TSLPSGDTFSSVGTHRFVQKPYSRSSMSIDLVSASDDVHRFSSQVEEMVQNHMNYSLQDVGGDANWQL
 VVEEGEMKVYRREVEENGIVLDPLKATHAVKGVGTGHEVCNYFWNVDRNDWETTIEFHVVELADNAII
 VYQTHKRVPASQRDVLVLSAIRKIPALTENDPETWIVCNFSVDHDSAPLNNRCVRAKINIAMICQTLVLS
 PPEGDQEISRDNILCKITYVANVNPGGWAPASVLRVAVAKREYPKFLKRFTSYVQEKTAGKPILF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_023420

ORF Size: 1872 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_023420.2 , NP_075909.1
RefSeq Size:	5363 bp
RefSeq ORF:	1875 bp
Locus ID:	68018
UniProt ID:	Q9EQG9
Cytogenetics:	13 D1
MW:	71.6 kDa
Gene Summary:	Shelters ceramides and diacylglycerol lipids inside its START domain and mediates the intracellular trafficking of ceramides and diacylglycerol lipids in a non-vesicular manner. [UniProtKB/Swiss-Prot Function]