

Product datasheet for **MR217782L4V**

Zc3h14 (NM_001008506) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zc3h14 (NM_001008506) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zc3h14
Synonyms:	1010001P15Rik; 1700016A15Rik; 2700069A02Rik; AU014748
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008506
ORF Size:	1812 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217782).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008506.2 , NP_001008506.2
RefSeq Size:	3163 bp
RefSeq ORF:	1815 bp
Locus ID:	75553
UniProt ID:	Q8BJ05
Cytogenetics:	12 E
Gene Summary:	Involved in poly(A) tail length control in neuronal cells. Binds the polyadenosine RNA oligonucleotides.[UniProtKB/Swiss-Prot Function]



[View online »](#)