

Product datasheet for **MR217783L3V**

Aldh3a1 (NM_001112725) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aldh3a1 (NM_001112725) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aldh3a1
Synonyms:	Ahd-4; Ahd4; Aldh; Aldh3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001112725
ORF Size:	1359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217783).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001112725.1 , NP_001106196.1
RefSeq Size:	1597 bp
RefSeq ORF:	1362 bp
Locus ID:	11670
UniProt ID:	P47739
Cytogenetics:	11 37.96 cM



[View online »](#)

Gene Summary:

ALDHs play a major role in the detoxification of alcohol-derived acetaldehyde (Probable). They are involved in the metabolism of corticosteroids, biogenic amines, neurotransmitters, and lipid peroxidation (Probable). Oxidizes medium and long chain aldehydes into non-toxic fatty acids (PubMed:25286108). Preferentially oxidizes aromatic aldehyde substrates (PubMed:11784860). Comprises about 50 percent of corneal epithelial soluble proteins (PubMed:11784860). May play a role in preventing corneal damage caused by ultraviolet light (PubMed:10376761).[UniProtKB/Swiss-Prot Function]