

## Product datasheet for **MR217795L4V**

### **Pdpx (NM\_020271) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pdpx (NM_020271) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pdpx
Synonyms:	1600027H05Rik; AB041662; PLPP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020271
ORF Size:	876 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217795).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020271.3</a> , <a href="#">NP_064667.2</a>
RefSeq Size:	2018 bp
RefSeq ORF:	879 bp
Locus ID:	57028
UniProt ID:	<a href="#">P60487</a>
Cytogenetics:	15 E1



[View online »](#)

**Gene Summary:**

Protein serine phosphatase that dephosphorylates 'Ser-3' in cofilin and probably also dephosphorylates phospho-serine residues in DSTN. Regulates cofilin-dependent actin cytoskeleton reorganization. Required for normal progress through mitosis and normal cytokinesis. Does not dephosphorylate phospho-threonines in LIMK1. Does not dephosphorylate peptides containing phospho-tyrosine (By similarity). Pyridoxal phosphate (PLP) phosphatase, which also catalyzes the dephosphorylation of pyridoxine 5'-phosphate (PNP) and pyridoxamine 5'-phosphate (PMP), with order of substrate preference PLP > PNP > PMP (By similarity) (PubMed:24338473, PubMed:24338687).[UniProtKB/Swiss-Prot Function]