

Product datasheet for **MR217822L3V**

Cd300lf (NM_001169153) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd300lf (NM_001169153) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd300lf
Synonyms:	CLIM1; CLM-1; CLM1; Digr2; F730004D16Rik; IgSF13; IREM1; LMIR3; Pigr3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_001169153
ORF Size:	1011 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217822).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001169153.1 , NP_001162624.1
RefSeq Size:	1885 bp
RefSeq ORF:	1014 bp
Locus ID:	246746
UniProt ID:	Q6SJQ7
Cytogenetics:	11 E2



[View online »](#)

Gene Summary:

Acts as an inhibitory receptor for myeloid cells and mast cells (PubMed:17438331). Positively regulates the phagocytosis of apoptotic cells (efferocytosis) via phosphatidylserine (PS) recognition; recognizes and binds PS as a ligand which is expressed on the surface of apoptotic cells (PubMed:21865548). Plays an important role in the maintenance of immune homeostasis, by promoting macrophage-mediated efferocytosis and by inhibiting dendritic cell-mediated efferocytosis (PubMed:26768664). Negatively regulates Fc epsilon receptor-dependent mast cell activation and allergic responses via binding to ceramide which acts as a ligand (PubMed:23123064). May act as a coreceptor for interleukin 4 (IL-4). Associates with and regulates IL-4 receptor alpha-mediated responses by augmenting IL-4- and IL-13-induced signaling (PubMed:26124135). Negatively regulates the Toll-like receptor (TLR) signaling mediated by MYD88 and TRIF through activation of PTPN6/SHP-1 and PTPN11/SHP-2 (By similarity). Inhibits osteoclast formation (PubMed:14662855). Induces macrophage cell death upon engagement (PubMed:18097021).[UniProtKB/Swiss-Prot Function]