

Product datasheet for **MR217831L4V**

Crlf3 (NM_018776) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Crlf3 (NM_018776) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Crlf3
Synonyms:	BB164954; Creme; Creme9; cyto; Cytor4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018776
ORF Size:	1326 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217831).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018776.2 , NP_061246.1
RefSeq Size:	2396 bp
RefSeq ORF:	1329 bp
Locus ID:	54394
UniProt ID:	Q9Z2L7
Cytogenetics:	11 47.43 cM



[View online »](#)

Gene Summary:

This gene encodes a cytokine receptor-like factor that contains a fibronectin type 3 domain. The encoded protein may act as a negative regulator of the cell cycle. There is a pseudogene for this gene on chromosome 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013]