

Product datasheet for **MR217833L3V**

Clasp1 (NM_177548) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Clasp1 (NM_177548) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Clasp1
Synonyms:	1700030C23Rik; 5730583A19Rik; B130045P17Rik; mKIAA0622
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177548
ORF Size:	4404 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217833).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177548.2 , NP_808216.2
RefSeq Size:	7623 bp
RefSeq ORF:	4407 bp
Locus ID:	76707
Cytogenetics:	1 E2.3



[View online »](#)

Gene Summary:

Microtubule plus-end tracking protein that promotes the stabilization of dynamic microtubules. Involved in the nucleation of noncentrosomal microtubules originating from the trans-Golgi network (TGN). Required for the polarization of the cytoplasmic microtubule arrays in migrating cells towards the leading edge of the cell. May act at the cell cortex to enhance the frequency of rescue of depolymerizing microtubules by attaching their plus-ends to cortical platforms composed of ERC1 and PHLDB2. This cortical microtubule stabilizing activity is regulated at least in part by phosphatidylinositol 3-kinase signaling. Also performs a similar stabilizing function at the kinetochore which is essential for the bipolar alignment of chromosomes on the mitotic spindle (By similarity).[UniProtKB/Swiss-Prot Function]